

# 360° Sanitizing & Disinfection System

Effectively Reducing Hospital-Acquired Infections.



- **3x more** surface coverage than traditional pump sprayers, rags or wipes
- **Faster, safer** and more consistent application
- **Touchless** electrostatic delivery of PURTABS
- Advanced battery technology – **4 hours continuous run**
- Portable and Cordless
- Easy for EVS and housekeeping to use

- EPA registered sanitizer & disinfectant tablets that readily dissolve in water; **Kills C. diff spores** in 4 minutes
- Cost effective and sustainable
- Eliminates and controls odors
- **Safer** for humans and materials
- **Stable** once diluted, and keeps potency and effectiveness for 7 days



Hard to reach places and soft surfaces



All touchpoints covered with electrostatic wrap

# PURTABS™

EFFERVESCENT SANITIZING / DISINFECTION TABLETS



## - EPA registered sanitizer & disinfectant multipurpose tablets that readily dissolve in water:

- NSF D2 Food Contact Surface Safe Sanitizer (100ppm)
- Sporicidal disinfectant; **Kills C. diff spores** in 4 minutes (4306ppm)

## - FAC85: Environmentally Preferable Cleaning Products

## - NIOSH/OSHA compliant with all safety regulations for exposure

## - A safer alternative to bleach in any situation. With zero compromises on cost, efficacy, or convenience

Kill Claims	Concentration	Dilution
Food contact surface sanitizer (No Rinse) 60 seconds contact time	100 ppm	(1) 334mg tab / 1qt (1) 3.3g tab / 2.5 gal (1) 13.1 tab / 10 gal
Hospital grade disinfectant claims (Bleach alt.) Cold & flu virus* (including H1N2) Salmonella & Staph.	538 ppm	(6) 334mg tab / 1 qt (1) 3.3g tab / 2 qts
Bloodborne pathogen, Herpes & Polio, Norovirus, MRSA & Ecoli, HIV & Hep B claims Many Animal pathogens	1076 ppm	(11) 334mg tab / 1 qt (1) 3.3g tab / 1 qt (1) 13.1 tab / 1 gal
C. Diff spore 10 minute claim	2153 ppm	(2) 3.3g tab / 1 qt (2) 13.1g tab / 1 gal
C. Diff spore 4 minute claim	4306 ppm	(4) 3.3 g tab / 1 qt (4) 13.1 g tab / 1 gal

# PROTEXUS™

ELECTROSTATIC SPRAYING SYSTEMS



## - Touchless disinfection

- In real-world use, Protexus Electrostatic Sprayers can coat all surfaces of a 500 sq/ft area with PURTABS in under 3 minutes.

- Greater than 80% time savings in any chemical application, disinfection or sanitizing process

### LIGHTWEIGHT & PORTABLE

Light enough at 3.8 lbs. to quickly navigate from surface and room to room

2 MINUTES

PROTEXUS ELECTROSTATIC SPRAYER

### QUIET

Being only 60 db., its great for environments where noise is important

### ADJUSTABLE NOZZLE

Allows the versatility of adjusting particle size for the required application

12 MINUTES

HAND PUMP TRIGGER SPRAYER

### VERSATILE

Ideal for sanitizing, disinfecting, odor control



Concentrated PURTABS tablets weigh **95% less** per case than liquid concentrate sanitizer/disinfectant. Making it **more ergonomic** to lift, and uses less plastics and cardboard, creating less solid waste.



PURTABS **saves money** on shipping, storage space and overall chemical cost.



PURTABS **weighs 7x less** than traditional liquid concentrate sanitizers and disinfectants, so **you pay 7x less** in shipping costs. PURTABS takes up less square footage, saving money on storage.



PURTABS is **less-toxic** and **less-corrosive** to humans, materials and the environment.

Mfr No.	Description	Qty
PX200ES	PROTEXUS CORDLESS ELECTROSTATIC SPRAYER	EA
ESPT3.3G	3.3G PURTABS - Hospital grade disinfectant (200 tabs/tub)	6 TUBS/CT
ESPT13.1G	13.1G PURTABS - Hospital grade disinfectant (256 tabs/tub)	2 TUBS/CT
ESPT334MG	334MG PURTABS - Food contact surface safe sanitizer (100 tabs/box)	12 BOXES/CT
PX30ES	REPLACEMENT TANK	EA



**HMIS**  
1 HEALTH  
0 FLAMMABILITY  
0 REACTIVITY