



# Proactive Sanitizing & Disinfection

Effectively Preventing Colds, Flu & Infections in Athletic Facilities



- **3x more** surface coverage
- **Fastest & easiest** way to disinfect facilities with more consistent application
- **Touchless** electrostatic delivery of PURTABS
- Advanced battery technology – **4 hours continuous run time**
- Portable and Cordless
- Easy to use & train staff

- EPA registered sanitizer & disinfectant tablets that readily dissolve in water
- **Safer** for humans, materials and the environment
- Eliminates and controls odors
- **Stable** once diluted, and keeps potency and effectiveness for 7 days
- Cost effective and sustainable
-  **D2 Certified Food Contact Surface Safe Sanitizer at 100 ppm**



# PURTABS™

EFFERVESCENT SANITIZING / DISINFECTION TABLETS



- EPA registered sanitizer & disinfectant multipurpose tablets that readily dissolve in water:
  - NSF D2 Food Contact Surface Safe Sanitizer (100ppm)
  - Kills Norovirus, MRSA, Cold and Flu viruses
- FAC85: Environmentally Preferable Cleaning Products
- NIOSH/OSHA compliant with all safety regulations for exposure
- A safer alternative to bleach in any situation. With zero compromises on cost, efficacy, or convenience


Kill Claims	Concentration	Dilution
Food contact surface sanitizer (No Rinse) 60 seconds contact time	100 ppm	(1) 334mg tab / 1qt (1) 3.3g tab / 2.5 gal (1) 13.1 tab / 10 gal
Hospital grade disinfectant claims (Bleach alt.) Cold & flu virus* (including H1N2) Salmonella & Staph.	538 ppm	(6) 334mg tab / 1 qt (1) 3.3g tab / 2 qts
Bloodborne pathogen, Herpes & Polio, Norovirus, MRSA & Ecoli, HIV & Hep B claims Many Animal pathogens	1076 ppm	(11) 334mg tab / 1 qt (1) 3.3g tab / 1 qt (1) 13.1 tab / 1 gal
C. Diff spore 10 minute claim	2153 ppm	(2) 3.3g tab / 1 qt (2) 13.1g tab / 1 gal
C. Diff spore 4 minute claim	4306 ppm	(4) 3.3 g tab / 1 qt (4) 13.1 g tab / 1 gal

# PROTEXUS™

ELECTROSTATIC SPRAYING SYSTEMS



- Touchless disinfection
- In real-world use, Protexus Electrostatic Sprayers can coat all surfaces of a 500 sq/ft area with PURTABS in under 3 minutes.
- Greater than 80% time savings in any chemical application, disinfection or sanitizing process



Concentrated PURTABS tablets weigh **95% less** per case than liquid concentrate sanitizer/disinfectant. Making it **more ergonomic** to lift, and uses less plastics and cardboard, creating less solid waste.



PURTABS **saves money** on shipping, storage space and overall chemical cost.



PURTABS **weighs 7x less** than traditional liquid concentrate sanitizers and disinfectants, so **you pay 7x less** in shipping costs. PURTABS takes up less square footage, saving money on storage.



PURTABS is **less-toxic** and **less-corrosive** to humans, materials and the environment.

Mfr No.	Description	Qty
PX200ES	PROTEXUS CORDLESS ELECTROSTATIC SPRAYER	EA
ESPT3.3G	3.3G PURTABS - Hospital grade disinfectant (200 tabs/tub)	6 TUBS/CT
ESPT13.1G	13.1G PURTABS - Hospital grade disinfectant (256 tabs/tub)	2 TUBS/CT
ESPT334MG	334MG PURTABS - Food contact surface safe sanitizer (100 tabs/box)	12 BOXES/CT
PX30ES	REPLACEMENT TANK	EA



**HMS**  
1 HEALTH  
0 FLAMMABILITY  
0 REACTIVITY